Moreton Group Medical Services

Comparison of Defibrillators





| | Mindray D1 Public | Philips Heart Start |
|--|----------------------|---------------------|
| 7-inch video display | \checkmark | 3C |
| A help button with video instruction guide | \checkmark | x |
| Bold and simple graphics help users focus on one task at a time | \checkmark | × |
| Time to shock from the end of CPR interval | Less than 10 seconds | Typically 8 seconds |
| Automatic volume adjustment and screen brightness adjustment | \checkmark | × |
| Optional dual language capability | \checkmark | × |
| High IP Rating for protection against water and dust | ✓ IP 55 | ≭ IP21 |
| Child mode button for quick access to reduced energy | \checkmark | × |
| AED to automatically detect a child without the use of child pads | \checkmark | × |
| Escalating biphasic full energy for difficult-to-defibrillate patients | 200– 360J | Only 150J fixed |
| Certified for use in Helicopter via Vibration Tests | \checkmark | × |
| Drop tested to 1.5m on 6 different angles | \checkmark | sc |
| Voice Recording in an event - programmable | \checkmark | ગ્ર |
| High-quality CPR | | |
| Metronome and CPR coaching to audibly guide userS | \checkmark | \checkmark |
| Visual CPR on screen feedback and instructions | \checkmark | ગ્ર |
| Programmable to most Emergency Protocols up to 360J | \checkmark | ઝર |
| CPR feedback, with the option of an accelerometer | \checkmark | ગ્ર |
| Connected. Ready. (Options) | | |
| AED program management system via Wi-Fi or cellular network | \checkmark | ગ |
| Location tracking of AED via a network | \checkmark | ઝર |
| Wireless monitoring of battery charge and self-tests | \checkmark | JC |
| Email notification via Wi-Fi of device in use | \checkmark | JC |
| Live ECG through the electrode pads onto the large 7 inch screen | \checkmark | ઝર |
| Live ECG event data for arriving emergency responders to view | \checkmark | ગ્ર |
| | | |

03 6278 3029

moretongroup.com